



TX



THE TILITE TX. STRONGER, LIGHTER AND DESIGNED TO OUT PERFORM!

WHY TITANIUM?

Reason No. 1: Titanium Absorbs Vibration.

Reason No. 2: Titanium Won't Break.

Reason No. 3: Titanium has the Highest Strength-to-Weight Ratio of Any Metal on Earth.

TILITE INNOVATION.

- Exclusive **Rigid Folding Mechanism "RFM"** increases performance by reducing flex in the side frames.
- **NEW Easy Fold RFM** makes folding the TX a breeze.
- **Bullet** caster housing reduces weight and gives the TX a stylish look.
- **Performance 5-Spoke Soft Roll Wheel** weighs just 1.6 ounces, 40% lighter than the original 3-Spoke Wheel.

TECHNICAL SPECIFICATIONS.

Frame Type	Titanium Folding – Minimal Adjustability
Frame Style	Fixed Front End
Total Width Folded*	10.75"
Total Chair Width*	23"
Weight Limit	250 lbs.
Seat Widths	12" - 20"
Seat Depths	14" - 20"
Front Seat Heights	15" - 21"
Rear Seat Heights	14" - 19"
Cambers	0°, 2° or 4°
Front Frame Angles	70°, 80° or 85°
Front Wheels	22 Options
Rear Wheels	61 Options
Handrims	12 Options
Center of Gravity Adjustment	Fixed or 3"
Transport Weight*	12.3 lbs.
Convenience	Folding Angle Adjustable Seat Back
Aesthetics	6 Titanium Frame Finishes 19 Frame Colors 4 Color Anodize Packages 7 Upholstery Colors 8 Spinerigy Wheel Colors
HCPCS Code	K0009
Warranty	Lifetime on Frame; Limited Lifetime on Components

* Weights and widths based on 16" x 16" chair with 0° camber, minimum rear wheel spacing and all standard features. Rear wheels are not included in transport weight. For reference, a pair of standard 24" rear wheels weighs 7.6 lbs.

TIFIT IS EXCLUSIVELY OFFERED BY TILITE.

The difference between a TiFit chair and all other chairs is that a TiFit chair is tailored specifically to you: Every frame is made to measure. No tube is cut, bent, or welded until our skilled fabricators have a full-scale drawing of your chair based on your exact measurements. Only TiLite offers TiFit. The optimal fit of a TiFit chair brings unparalleled performance and function. You will feel the difference in everything you do.

CUSTOM CONFIGURATION THE TILITE WAY.

Custom configuration is when you are able to select options and components on a chair to meet your needs. While other companies offer some custom configuration, TiLite offers the widest range of options from which to choose. More choices mean more ability to truly create a custom configuration of options that best meets your needs. Choice matters and that is why TiLite offers more options than any other manufacturer! The choice is yours.



TiLite | 800.545.2266 | tilite.com
© 2011 TiSport, LLC.

All TiLite wheel chairs are shown with available options. CF_TX_0911RevC